

### (g) Simultaneous Measurement of NADPH Oxidase Activity and Phagosomal Oxidation with Hydroethidine and 2',7'-Dichlorofluorescein Diacetate

**Submitted by:**

Gregor Rothe and Günter Valet  
Arbeitsgruppe Zellbiochemie  
Max-Planck-Institut für Biochemie  
W-8033 Martinsried  
Germany

Phone: 49-(89)-8578-2525  
Fax: 49-(89)-8578-3777

**Outline**

*During the oxidative burst, phagocytic cells release superoxide anion through the membrane-bound NADPH oxidase. Hydrogen peroxide, produced by dismutation of superoxide anion, is the substrate for the myeloperoxidase-catalyzed oxidation inside the phagosome. Hydroethidine can be oxidized to red fluorescent ethidium bromide by superoxide anion directly, whereas the oxidation of 2',7'-dichlorofluorescein to green fluorescent 2',7'-dichlorofluorescein requires hydrogen peroxide and peroxidase. The simultaneous two-color analysis of hydroethidine oxidation and 2',7'-dichlorofluorescein oxidation permits the detection of selective defects of the phagosomal oxidation following lysosomal degranulation in diseases such as sepsis.*

**Specimen:** 3 ml heparinized human blood (10 U heparin/ml)

**Reagents**

HBSS without phenol red or bicarbonate, supplemented with 10 mM HEPES (pH 7.35)

hydroethidine (MW 315)

stock solution: 10 mM in DMF (3.15 mg/ml)

working solution: 1 mM (1:10 dilution of stock in HBSS)

2',7'-dichlorofluorescein diacetate (DCFH-DA) (MW 487.3)

stock solution: 10 mM (4.87 mg/ml) in DMF

working solution: 1 mM (1:10 dilution of stock in HBSS)

propidium iodide (PI) (MW 668.4)

stock solution: 3 mM (2 mg/ml) in 5 mM HEPES-buffered saline (0.15 M NaCl, pH 7.35)

phorbol 12-myristate 13-acetate (PMA) (MW 616.8)

stock solution: 1 mM in DMF

working solution: 10  $\mu$ M (1:100 dilution of stock in HBSS)

**Procedure**

1. Layer 3 ml heparinized blood carefully on top of 3 ml lymphocyte separation medium. Allow erythrocytes to sediment for 40 minutes at room temperature.
2. Withdraw the upper 800  $\mu$ l supernatant plasma and store on ice. This will contain platelets and approximately  $2 \times 10^7$ /ml unseparated leukocytes.
3. Put 1.00 ml HBSS, 20  $\mu$ l cell suspension, and 10  $\mu$ l hydroethidine working solution in a 12 x 75 mm polypropylene test tube (final hydroethidine concentration 10  $\mu$ M). Incubate for 5 minutes at 37°C. Add 10  $\mu$ l PMA working solution (final PMA concentration 100 nM). Continue incubation, taking 250  $\mu$ l aliquots at 10, 20, and 30 minutes after addition.
4. Counterstain dead cells by incubating 250  $\mu$ l stained cell suspension with 5  $\mu$ l 3mM PI for 3 minutes on ice (final PI concentration 60  $\mu$ M).
5. Run on flow cytometer.

Excitation: 488 nm (argon laser) or high pressure mercury arc lamp with 470-500 nm bandpass filter

Filters: 510-530 bandpass for 2',7'-dichlorofluorescein green fluorescence

600 nm longpass for ethidium bromide (viable cells) and PI (dead cells) red fluorescence



# HANDBOOK

*of*

# FLOW

*Edited by J. Paul Robinson*  
*Director, Purdue University*  
*Cytometry Laboratories*  
*Hansen Life Sciences Research Building*  
*Purdue University*  
*West Lafayette, Indiana*

# CYTOMETRY

# METHODS

## *Associate Editors*

*Zbigniew Darzynkiewicz*  
*New York Medical College*  
*Valhalla, New York*

*Phillip Dean*  
*Lawrence Livermore*  
*National Laboratories*  
*Livermore, California*

*Lynn Dressler*  
*University of Arizona*  
*Tucson, Arizona*

*Hans Tanke*  
*Sylvius Laboratorium*  
*Leiden, The Netherlands*

*Leon Wheeless*  
*University of Rochester*  
*Rochester, New York*

 **WILEY-LISS**

A JOHN WILEY & SONS, INC., PUBLICATION  
New York • Chichester • Brisbane • Toronto • Singapore

**Address all Inquiries to the Publisher**  
**Wiley-Liss, Inc., 605 Third Avenue, New York, NY 10158-0012**

**Copyright © 1993 Wiley-Liss, Inc.**

**Printed in the United States of America.**

Under the conditions stated below the owner of copyright for this book hereby grants permission to users to make photocopy reproductions of any part or all of its contents for personal or internal or organizational use, or for personal or internal use of specific clients. This consent is given on the condition that the copier pay the stated per-copy fee through the Copyright Clearance Center, Incorporated, 27 Congress Street, Salem, MA 01970, as listed in the most current issue of "Permissions to Photocopy" (Publisher's Fee List, distributed by CCC, Inc.), for copying beyond that permitted by sections 107 or 108 of the US Copyright Law. This consent does not extend to other kinds of copying, such as copying for general distribution, for advertising or promotional purposes, for creating new collective works, or for resale.

While the authors, editors, and publisher believe that drug selection and dosage and the specification and usage of equipment and devices, as set forth in this book, are in accord with current recommendations and practice at the time of publication, they accept no legal responsibility for any errors or omissions, and make no warranty, express or implied, with respect to material contained herein. In view of ongoing research, equipment modifications, changes in governmental regulations and the constant flow of information relating to drug therapy, drug reactions, and the use of equipment and devices, the reader is urged to review and evaluate the information provided in the package insert or instructions for each drug, piece of equipment, or device for, among other things, any changes in the instructions or indication of dosage or usage and for added warnings and precautions.

The publication of this volume was facilitated by the authors and editors, who submitted the text in a form suitable for direct reproduction without subsequent editing or proofreading by the publisher.

#### **Library of Congress Cataloging-in-Publication Data**

Handbook of flow cytometry methods / editor, J. Paul Robinson ;  
associate editors, Zbigniew Darzynkiewicz ... [et al.].

p. cm.

Includes bibliographical references (p. ) and index.

ISBN 0-471-59634-5

1. Flow cytometry—Laboratory manuals. I. Robinson, J. Paul.

QH585.5.F56H36 1993

574.87'028—dc20

92-47082

CIP

**The text of this book is printed on acid-free paper.**